



MULTI-PARAMETER PATIENT SIMULATOR



PS-2105

USER MANUAL

**BC BIOMEDICAL
PS-2105
TABLE OF CONTENTS**

WARNINGS, CAUTIONS, NOTICES	ii
DESCRIPTION.....	1
LAYOUT.....	4
ECG--NORMAL SINUS RHYTHM.....	11
ECG--ARRHYTHMIAS	13
ECG--PACEMAKER.....	20
ECG--PERFORMANCE	22
TRAINING	24
SETUP	29
POWER UP SETTINGS.....	31
SpO ₂ (OPTION).....	33
OUTPUT CONNECTIONS	34
MANUAL REVISIONS.....	35
LIMITED WARRANTY.....	35
SPECIFICATIONS	36
NOTES.....	38

This User Manual covers the following units:

- PS-2105

WARNING - USERS

The PS-2105 is for use by skilled technical personnel only.

WARNING - USE

The PS-2105 is intended for testing only and should never be used in diagnostics, treatment or any other capacity where it would come in contact with a patient.

WARNING - CONNECTIONS

All connections to patients must be removed before connecting the DUT to the PS-2105. A serious hazard may occur if the patient is connected when testing with the PS-2105.

CAUTION - MODIFICATIONS

The PS-2105 is intended for use within the published specifications. Any application beyond these specifications or any unauthorized user modifications may result in hazards or improper operation.

CAUTION - SERVICE

The PS-2105 is intended to be serviced only by authorized service personnel. Troubleshooting and service procedures should only be performed by qualified technical personnel.

CAUTION - INSPECTION

The PS-2105 should be inspected before each use for obvious signs of abuse or wear. The PS-2105 should not be used and should be serviced if any parts are in question.

CAUTION - CLEANING

Do not immerse. The PS-2105 should be cleaned by wiping gently with a damp, lint-free cloth. A mild detergent can be used if desired.

CAUTION - LIQUIDS

Do not submerge or spill liquids on the PS-2105. Do not operate the PS-2105 if it may have been exposed to fluid.

CAUTION - ENVIRONMENT

Exposure to environmental conditions outside the specifications can adversely affect the performance of the PS-2105. Allow the PS-2105 to acclimate to specified conditions for at least 30 minutes before attempting to operate it.



NOTICE – CE



The PS-2000 Series Simulators bear the  mark
Based on the following testing standards:

ELECTROMAGNETIC COMPATIBILITY DIRECTIVE

**EMC – Directive 89/336/EEC as amended by 92/31/EEC and 93/68/EEC
& Directive 91/263/EEC[TTE/SES]**

**EN 61326-1:1997 + A1:1998 + A2:2001 + A3:2003
“Electrical equipment for measurement, control and
laboratory use – EMC requirements”**

This equipment has been type tested by an independent, accredited testing laboratory
and compliance was demonstrated to the above standard to the extent applicable.

EMISSIONS
Radiated Emissions

EN 61326:1997 Annex C

IMMUNITY– CLASS C

EN 61000-4-2:1995
EN 61000-4-3:2006

Electrostatic Discharge
Radiated Electric Field Immunity

LOW VOLTAGE DIRECTIVE **EC – Directive 73/23/EC**

**EN 61010-1:2001
“Safety requirements for electrical equipment for measurement, control, and
laboratory use – General requirements”**

This equipment has been type tested and compliance was demonstrated
to the above standard to the extent applicable.

NOTICE – SYMBOLS

<u>Symbol</u>	<u>Description</u>
---------------	--------------------



Caution
(Consult Manual for Further Information)



Center Negative



Per European Council Directive 2002/95/EC,
do not dispose of this product as unsorted
municipal waste.

NOTICE – ABBREVIATIONS

AHA	American Heart Association
ANSI	American National Standards Institute
BPM	Beats Per Minute
C	Celsius
°	degree(s)
DUT	Device Under Test
ECG	Electrocardiogram
F	Fahrenheit
Hz	hertz
IEC	International Electrotechnical Commission
Lbs	pounds
LED	Light Emitting Diode
mm	millimeter(s)
mV	millivolt(s)
NEDA	National Electronic Distributors Association
USA	United States of America
VDC	Volts Direct Current

NOTICE – DISCLAIMER

BC GROUP INTERNATIONAL, INC. WILL NOT BE RESPONSIBLE FOR ANY INJURIES SUSTAINED DUE TO UNAUTHORIZED EQUIPMENT MODIFICATIONS OR APPLICATION OF EQUIPMENT OUTSIDE OF THE PUBLISHED INTENDED USE AND SPECIFICATIONS.

NOTICE – DISCLAIMER

BC GROUP INTERNATIONAL, INC. RESERVES THE RIGHT TO MAKE CHANGES TO ITS PRODUCTS OR SPECIFICATIONS AT ANY TIME, WITHOUT NOTICE, IN ORDER TO IMPROVE THE DESIGN OR PERFORMANCE AND TO SUPPLY THE BEST POSSIBLE PRODUCT. THE INFORMATION IN THIS MANUAL HAS BEEN CAREFULLY CHECKED AND IS BELIEVED TO BE ACCURATE. HOWEVER, NO RESPONSIBILITY IS ASSUMED FOR INACCURACIES.

NOTICE – CONTACT INFORMATION

BC BIOMEDICAL
BC GROUP INTERNATIONAL, INC.
3081 ELM POINT INDUSTRIAL DRIVE
ST. CHARLES, MO 63301
USA

1-800-242-8428

1-314-638-3800

www.bcgrouptl.com

sales@bcgrouptl.com

BC BIOMEDICAL PS-2105 PATIENT SIMULATOR

The Model PS-2105 is a Microprocessor based Patient Simulator. It provides ECG Simulation with 36 arrhythmias, 17 waveforms with constant QRS duration, 15 machine performance-testing waveforms and 7 paced rhythms, plus a unique training mode and an optional SpO₂ Output.

The PS-2105 makes viewing and selecting the desired waveforms and parameters quick and intuitive, with all operational information being available at the same time on a cursor-based graphic display, allowing for easy maneuvering through parameters and scrolling through available options.

The following are highlights of some of the main features:

PS-2105 (BASIC FEATURES):

- SIMPLE TO OPERATE
- NO CODES TO REMEMBER OR ENTER
- GRAPHICS DISPLAY WITH SIMULTANEOUS DETAILED STATUS OF PARAMETERS AND SCROLLING CONTROL OF OPTIONS
- DROP DOWN CHOICE SCREENS LIST ALL OPTIONS FOR PARAMETERS
- SPECIAL POWER UP FEATURE ALLOWS THE USER TO CHOOSE TO USE DEFAULT, LAST OR CUSTOM SETTINGS
- AUTO SEQUENCES FOR BPM AND PERFORMANCE
- 10 UNIVERSAL PATIENT LEAD CONNECTORS
- 9 VOLT BATTERY POWER
- % BATTERY LIFE DISPLAY
- LOW BATTERY INDICATOR
- DISPLAY BACKLIGHT
- FLASH PROGRAMMABLE FOR UPGRADES

ECG FUNCTIONS

The unit can produce a wide variety of ECG simulations. The user simply selects the parameters that match the desired output. Available settings:

- RATE: 30, 40, 45, 60, 80, 90, 100, 120, 140, 160, 180, 200, 220, 240, 260, 280, 300 BPM
- AMPLITUDE:
0.5, 1.0, 1.5, 2.0 mV (Lead II)
- S-T SEGMENT ELEVATION:
± 0, 0.05, 0.1, 0.2, 0.3, 0.4, 0.5, 0.6, 0.7, 0.8 mV.
- ARTIFACTS:
50 HZ, 60 HZ, MUSCLE, BASELINE WANDER, RESPIRATION

- QRS INTERVAL:
ADULT (80 ms) OR PEDIATRIC (40 ms)
- AUTOMATIC MODE

ECG-PERFORMANCE FUNCTIONS

The unit will generate Sine, Square, Triangular and Pulse waveforms with adjustable amplitudes for performance testing. A special Automatic mode is available to auto sequence through the entire range of waveforms. Available settings:

- SINE:
0.1, 0.5, 5, 10, 40, 50, 60, 100 Hz
- SQUARE:
0.125, 2 Hz
- TRIANGLE:
2, 2.5 Hz
- PULSE:
30, 60, 120 BPM; 60 ms WIDTH
- AMPLITUDE:
0.5, 1.0, 1.5, 2.0 mV (Lead II)
- AUTOMATIC MODE

PACEMAKER FUNCTIONS

Seven different pacemaker waveforms may be simulated. Additionally, the width and amplitude of the pacer pulse may be selected. Available settings:

- WAVEFORMS:
ATRIAL PACER, ASYNCHRONOUS, NON-CAPTURE,
NON-FUNCTION, DEMAND - OCCASIONAL,
DEMAND - FREQUENT, AV - SEQUENTIAL
- PULSE HEIGHT:
1, 2, 5, 10 mV
- PULSE WIDTH:
0.1, 0.5, 1.0, 1.5, 2.0 ms

ARRHYTHMIA FUNCTIONS

The unit can simulate 36 different arrhythmias. For ease of selection, they are arranged into four basic groups. Where applicable, both manual and automatic triggering of the waveform is available. Available settings:

- 36 DIFFERENT ARRHYTHMIAS
- FOUR GENERAL GROUPS:
SUPRAVENTRICULAR
PREMATURE
VENTRICULAR
CONDUCTION

MANUAL AND AUTOMATIC TRIGGERING

TRAINING

The unit has a special training mode that may be used to aid users in practicing the identification of arrhythmias. A series of settings allows the feature to be customized to fit the exact training requirement. Available settings:

- TIMER:
MANUAL, 10, 15, 20, 25, 30 SEC
- RANDOMIZER:
OFF, ON
- ARRHYTHMIAS:
ALL, SUBSET

SPO2 SIMULATION (Option)

When used with the MSP-2100 external module and FingerSim family of SpO₂ finger simulators, the system will provide a pulse synchronized SpO₂ output for NSR rates.

Available settings:

- RATE: 30, 40, 45, 60, 80, 90, 100, 120, 140, 160, 180 BPM
- SpO₂ OUTPUT:
80, 90, 97 %

LAYOUT

This section looks at the layout of a PS-2105 and gives descriptions of the elements that are present.

10 Universal Patient Lead Connectors:

RA	R
LA	L
RL(-)	N
LL	F
V1	C1
V2	C2
V3	C3
V4	C4
V5	C5
V6	C6

Graphical LCD Display:

Shows Parameters for ECG, including BPM, Adult/Pediatric, Amplitude, Artifact and ST Segment

7 Light Touch Keys for Selecting Parameters and Settings:

LEFT and RIGHT Curved Arrows for Moving through Parameters
 UP and DOWN Arrows for Scrolling through Options
 ENTER for Selecting Option
 CHOICES for Displaying Submenu of All Options for a Given Parameter
 QUIT for Returning to Previous Status

7-Pin Mini-DIN Plug Connector for Aux Functions

Battery Eliminator Receptacle


Power Key for Turning Unit On and Off

6 Light Touch Keys for Category Selection:

Normal Sinus Rhythm
 Arrhythmias
 Pacemaker
 Performance
 Training
 Setup







9V Battery Compartment (Back)





Back Light Key for turning on and off the backlight












General Operation

The unit is controlled by 15 light touch keys. They allow the user to move around within the displayed parameters, select the desired options, choose a specific category and control the setup and power for the unit. When a key is depressed there is an audio click when it is accepted, or a razz tone if the key is invalid.


A graphical LCD display provides the user with information about the current status of the ECG settings. The   keys move the block cursor through the displayed information; highlighting the parameter available for selection. The   keys change the options for the highlighted parameter. The cursor begins flashing if the parameter has been changed. The  key selects the changed option. The  key returns to the previous state, without any changes being made.


To make option selection even easier and to make memorizing and using codes unnecessary, the  key will bring up a screen that displays all the options for the selected parameter. The   and  keys can then be used to quickly scroll through the available options and select the desired setting.


Five category keys allow for quick setting of output waveforms. The   ,  ,  and  keys move the display directly to the selected category. The   or  keys can then be used to scroll through and select the desired settings.


The  key opens a screen that allows the user to select the unit's general output settings, as well as setup for the system.


Category Keys

The  key enters the NSR category.


The  key enters the arrhythmia category and changes the first line in the display to the first arrhythmia choice.

The  key enters the pacemaker category and changes the first line in the display to the first pacemaker waveform choice

The  key enters the machine performance-testing category and changes the first line in the display to the first performance waveform choice.

The  key opens a screen that allows the user to set the conditions for and start the special training mode.


Power Key

The  key turns the unit on and off. To turn off the unit, the key must be held for 1 second.

Backlight

The Graphic LCD display may be viewed with or without the backlight. Depressing any key will activate the backlight. However, since the backlight will drain the battery if left on, it will automatically shut off after a few seconds when running on battery power. (Note: This time is selectable in the System Setup screen).

The intensity of the backlight can be adjusted in the System Setup screen to conserve battery life.

The  key is provided to toggle the backlight on or off at any time.

NOTE: The backlight parameter in the System Setup screen may be set to Off, 1-30 sec Timed or Manual.

ECG Waveforms

The microprocessor has stored in its memory all of the digitalized waveforms. It sends the individual lead waveforms to D/A converters, which generate accurate analog representations. The waveforms are then sent through resistor networks, developing the appropriate signals on the output terminals.

Universal Patient Lead Connectors

The 10 Universal Patient Lead Connectors allow for 12-lead ECG simulation with independent outputs. AHA and IEC color-coded labels are located on the face of the unit to aid in connecting the corresponding U.S. and International Patient Leads.

AHA Label	IEC Label	Description
RA	R	Right Arm
LA	L	Left Arm
RL	N	Right Leg (reference or ground)
LL	F	Left Leg
V1 V2 V3 V4 V5 V6	C1 C2 C3 C4 C5 C6	V Leads (U.S. Canada), also referred to as pericardial, precordial or unipolar chest leads Chest Leads (International)


High Level Output (+)

A high level ECG output signal (200 x Amplitude Setting) is available in the Aux 7-Pin mini-DIN connector.

Auto Power Off

The unit may be programmed to automatically turn off after a selected number of minutes of no key activity to conserve the battery. (Note: This time is selectable in the System Setup screen).

Battery

The unit utilizes two 9 Volt Alkaline Batteries in the rear battery compartment. When the unit detects a LOW BATTERY condition (5% Battery Life), a warning window will appear once per minute to alert the user. The  key may be used to clear this window and continue use of the unit. If the battery is not replaced before the battery reaches a critical level (0 % Battery Life), the unit will shut down. (The percentage of life left in the batteries can be viewed in the System Setup screen.)

Battery Eliminator



The unit has a 2.1 mm jack for connecting a 9-Volt Battery Eliminator (Optional). Note: The Battery Eliminator will not charge the battery.

Power Up Settings

The unit may be setup to turn on using either the factory default settings, the same settings that it had when last turned off or a custom set of parameters as previously saved by the user (See Power Up Settings section for details).

Automatic Modes


The ECG NSR Rate, ECG Performance and Static Blood Pressure Parameters all allow for an automatic setting. In each of these, the unit will sequence through the full range of settings automatically at a fixed rate (as selected in the Auto Step Time Parameter). When in this mode, the time remaining in each step is displayed.

The  key may be used to manually advance to the next step. The  key is used to terminate the mode.


ECG – NORMAL SINUS RHYTHM






The PS-2105 can send waveforms to ECG machines in 3, 5 or 12-lead configurations. It has independent outputs for each signal lead, referenced to the right leg.





Normal Sinus Rhythm (NSR) occurs when the heartbeat is normal, beating at a rate between 50 and 100 BPM with a standard QRS waveform shape and height. The PS-2105 simulates the NSR with a default pulse of 80 BPM, amplitude of 1.0 mV on Lead II, P-R interval of 160 milliseconds, no Artifact and no ST Segment elevation.

The PS-2105 is placed into NSR mode by pressing the  category key.

The display will resemble the following:



The rate, amplitude, adult/pediatric, artifact and ST elevation or depression can be selected by using   to highlight the parameter to change and using   to scroll to the desired value. Then  is used to accept the new setting.


Alternately, to see a submenu of all the options for a highlighted parameter, use . Use   to scroll to the desired option. Then  is used to accept the new setting.

Auto Rate

If the BPM parameter is set to AUTO, the unit will automatically sequence through all of the BPM settings, starting with 30 BPM, incrementing at a fixed interval. The interval may be set in the System Setup Menu under “Auto Step Time”.

Displays time (seconds) remaining before advancing to next rate.


Normal Sinus Rhythm	
30 BPM(29)	Adult
Amplitude:	1.0 mV
Artifact:	none
ST Segment:	0.0 mV

The  key can be used to exit the Auto Mode during the sequence.

NOTE: ST Elevation or Depression is only active in Adult NSR at or below 180 BPM.






ECG – ARRHYTHMIAS





The PS-2105 can send arrhythmia waveforms to ECG machines in 3, 5 or 12-lead configurations. It has independent outputs for each signal lead, referenced to the right leg.


There are 36 Arrhythmias available that model abnormal heartbeats. The PS-2105 is placed into ARRHYTHMIA mode by pressing the  category key.

The display will resemble the following:


GROUP	ARRHYTHMIA																																																		
Premature	Atrial PAC - Auto	<div style="border: 1px solid black; padding: 5px; margin-bottom: 5px;"> <p>Premature Arrhythmias:</p> <p style="text-align: center;">Atrial PAC - Auto</p> <hr/> <p>Amplitude: 1.0 mV</p> <hr/> <p>Artifact: 60 Hz</p> </div> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 20%;">Ventricular</td> <td>Pair of PVCs - Auto</td> </tr> <tr> <td></td> <td>Pair of PVCs - Man</td> </tr> <tr> <td></td> <td>Run of 5 PVCs - Auto</td> </tr> <tr> <td></td> <td>Run of 5 PVCs - Man</td> </tr> <tr> <td></td> <td>Run of 11 PVCs - Auto</td> </tr> <tr> <td></td> <td>Run of 11 PVCs - Man</td> </tr> <tr> <td></td> <td>6 PVCs per Min</td> </tr> <tr> <td></td> <td>12 PVCs per Min</td> </tr> <tr> <td></td> <td>24 PVCs per Min</td> </tr> <tr> <td></td> <td>Freq Multifocal PVCs</td> </tr> <tr> <td></td> <td>Bigeminy</td> </tr> <tr> <td></td> <td>Trigeminy</td> </tr> <tr> <td></td> <td>Vent Tach</td> </tr> <tr> <td></td> <td>Vent Fib – Coarse</td> </tr> <tr> <td></td> <td>Vent Fib – Fine</td> </tr> <tr> <td></td> <td>Asystole</td> </tr> <tr> <td rowspan="9">Supraventricular</td> <td>Atrial Fib - Coarse</td> <td rowspan="9" style="vertical-align: top;"> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th colspan="2">AMPLITUDE</th> </tr> <tr> <td>0.5 mV</td> </tr> <tr> <td>1.0 mV*</td> </tr> <tr> <td>1.5 mV</td> </tr> <tr> <td>2.0 mV</td> </tr> </table> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th colspan="2">ARTIFACT</th> </tr> <tr> <td>None*</td> </tr> <tr> <td>50 Hz</td> </tr> <tr> <td>60 Hz</td> </tr> <tr> <td>Muscle</td> </tr> <tr> <td>Baseline Wander</td> </tr> <tr> <td>Respiration</td> </tr> </table> <p style="font-size: small; border: 1px dashed black; padding: 2px; margin-top: 5px;">* Indicates Default Setting (See Power Up Settings)</p> </td></tr></table>	Ventricular	Pair of PVCs - Auto		Pair of PVCs - Man		Run of 5 PVCs - Auto		Run of 5 PVCs - Man		Run of 11 PVCs - Auto		Run of 11 PVCs - Man		6 PVCs per Min		12 PVCs per Min		24 PVCs per Min		Freq Multifocal PVCs		Bigeminy		Trigeminy		Vent Tach		Vent Fib – Coarse		Vent Fib – Fine		Asystole	Supraventricular	Atrial Fib - Coarse	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th colspan="2">AMPLITUDE</th> </tr> <tr> <td>0.5 mV</td> </tr> <tr> <td>1.0 mV*</td> </tr> <tr> <td>1.5 mV</td> </tr> <tr> <td>2.0 mV</td> </tr> </table> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th colspan="2">ARTIFACT</th> </tr> <tr> <td>None*</td> </tr> <tr> <td>50 Hz</td> </tr> <tr> <td>60 Hz</td> </tr> <tr> <td>Muscle</td> </tr> <tr> <td>Baseline Wander</td> </tr> <tr> <td>Respiration</td> </tr> </table> <p style="font-size: small; border: 1px dashed black; padding: 2px; margin-top: 5px;">* Indicates Default Setting (See Power Up Settings)</p>	AMPLITUDE		0.5 mV	1.0 mV*	1.5 mV	2.0 mV	ARTIFACT		None*	50 Hz	60 Hz	Muscle	Baseline Wander	Respiration
	Ventricular		Pair of PVCs - Auto																																																
			Pair of PVCs - Man																																																
			Run of 5 PVCs - Auto																																																
			Run of 5 PVCs - Man																																																
			Run of 11 PVCs - Auto																																																
			Run of 11 PVCs - Man																																																
			6 PVCs per Min																																																
			12 PVCs per Min																																																
			24 PVCs per Min																																																
			Freq Multifocal PVCs																																																
			Bigeminy																																																
			Trigeminy																																																
			Vent Tach																																																
			Vent Fib – Coarse																																																
			Vent Fib – Fine																																																
			Asystole																																																
	Supraventricular		Atrial Fib - Coarse	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th colspan="2">AMPLITUDE</th> </tr> <tr> <td>0.5 mV</td> </tr> <tr> <td>1.0 mV*</td> </tr> <tr> <td>1.5 mV</td> </tr> <tr> <td>2.0 mV</td> </tr> </table> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th colspan="2">ARTIFACT</th> </tr> <tr> <td>None*</td> </tr> <tr> <td>50 Hz</td> </tr> <tr> <td>60 Hz</td> </tr> <tr> <td>Muscle</td> </tr> <tr> <td>Baseline Wander</td> </tr> <tr> <td>Respiration</td> </tr> </table> <p style="font-size: small; border: 1px dashed black; padding: 2px; margin-top: 5px;">* Indicates Default Setting (See Power Up Settings)</p>	AMPLITUDE		0.5 mV	1.0 mV*	1.5 mV	2.0 mV	ARTIFACT		None*	50 Hz	60 Hz	Muscle	Baseline Wander	Respiration																																	
AMPLITUDE																																																			
0.5 mV																																																			
1.0 mV*																																																			
1.5 mV																																																			
2.0 mV																																																			
ARTIFACT																																																			
None*																																																			
50 Hz																																																			
60 Hz																																																			
Muscle																																																			
Baseline Wander																																																			
Respiration																																																			

The grouping, arrhythmias and amplitude can be selected by using   to highlight the parameter to change and using   to scroll to the desired option. Then  is used to accept the new setting.

Alternately, to see a submenu of all the options for a highlighted parameter, use  Use   to scroll to the desired option. Then  is used to accept the new setting.

NOTE: While in the Arrhythmia Group choice screen, the  key may be used for a second time to jump directly to the arrhythmias choices for that group.

Auto/Manual

There are 12 arrhythmias that have both Automatic and Manual versions. Both versions output the same waveform; however, in the Manual version, the arrhythmia is triggered each time  is depressed. In the Auto versions, the arrhythmia is automatically triggered periodically.

The following is a brief description of how the PS-2105 simulates the available Arrhythmias:

PREMATURE		
Abbreviation	Arrhythmia	Description
Atrial PAC – Auto	Premature Atrial Contraction	NSR of 80 BPM with Periodic Abnormal 25 % early P waves (PAC, 7 NSR) (Continuous)
Atrial PAC – Man	Premature Atrial Contraction	NSR of 80 BPM with Periodic Abnormal 25 % early P waves (One-Time event)
Nodal PNC – Auto	Premature Nodal Contraction	NSR of 80 BPM with Periodic Abnormal 25 % early Nodal beat (PNC, 7 NSR) (Continuous)
Nodal PNC – Man	Premature Nodal Contraction	NSR of 80 BPM with Periodic Abnormal 25 % early Nodal beat (One-Time event)
PVC 1 – Auto	Standard Type 1 Premature Ventricular Contraction	NSR of 80 BPM with periodic left focus premature ventricular beats with 20% premature timing (PVC Type 1, 9 NSR) (Continuous)
PVC 1 – Man	Standard Type 1 Premature Ventricular Contraction	NSR of 80 BPM with periodic left focus premature ventricular beats with 20% premature timing (One-Time event)
PVC 1 Early - Auto	Early Type 1 Premature Ventricular Contraction	NSR of 80 BPM with periodic left focus premature ventricular beats with 33% premature timing (PVC Type 1, 9 NSR) (Continuous)
PVC 1 Early - Man	Early Type 1 Premature Ventricular Contraction	NSR of 80 BPM with periodic left focus premature ventricular beats with 33% premature timing (One-Time event)

PVC 1 R on T – Auto	R on T Type 1 Premature Ventricular Contraction	NSR of 80 BPM with periodic left focus premature ventricular beats with 65% premature timing, placing R on the previous T (PVC Type 1, 9 NSR) (Continuous)
PVC 1 R on T – Man	R on T Type 1 Premature Ventricular Contraction	NSR of 80 BPM with periodic left focus premature ventricular beats with 65% premature timing, placing R on the previous T (One-Time event)
PVC 2 – Auto	Standard Type 2 Premature Ventricular Contraction	NSR of 80 BPM with periodic right focus premature ventricular beats with 20% premature timing (PVC Type 2, 9 NSR) (Continuous)
PVC 2 – Man	Standard Type 2 Premature Ventricular Contraction	NSR of 80 BPM with periodic right focus premature ventricular beats with 20% premature timing (One-Time event)
PVC 2 Early - Auto	Early Type 2 Premature Ventricular Contraction	NSR of 80 BPM with periodic right focus premature ventricular beats with 33% premature timing (PVC Type 2, 9 NSR) (Continuous)
PVC 2 Early - Man	Early Type 2 Premature Ventricular Contraction	NSR of 80 BPM with periodic right focus premature ventricular beats with 33% premature timing (One-Time event)
PVC 2 R on T – Auto	R on T Type 2 Premature Ventricular Contraction	NSR of 80 BPM with periodic right focus premature ventricular beats with 65% premature timing, placing R on the previous T (PVC Type 2, 9 NSR) (Continuous)
PVC 2 R on T – Man	R on T Type 2 Premature Ventricular Contraction	NSR of 80 BPM with periodic right focus premature ventricular beats with 65% premature timing, placing R on the previous T (One-Time event)
Multifocal PVCS – Auto	Multifocal Premature Ventricular Contraction	NSR of 80 BPM with Type 1 and Type 2 PVCs (PVC Type 1, 2 NSR, PVC Type 2, 2 NSR) (Continuous)
Multifocal PVCS – Man	Multifocal Premature Ventricular Contractions	NSR of 80 BPM with Type 1 and Type 2 PVCs (PVC Type 1, 2 NSR, PVC Type 2) (One-Time event)

SUPRAVENTRICULAR		
Abbreviation	Arrhythmia	Description
Atrial Fib – Coarse	Atrial Fibrillation	Absence of P-wave, irregular P-R interval rate and a high level signal (Continuous)
Atrial Fib – Fine	Atrial Fibrillation	Absence of P-wave, irregular P-R interval rate and a low level signal (Continuous)
Atrial Flutter	Atrial Flutter	Repeating sequence of 5 atrial beats and 1 ventricular beat for twelve seconds, followed by a repeating sequence of 3 atrial beats and 1 ventricular beat for six seconds, followed by a repeating sequence of 2 atrial beats and 1 ventricular beat for six seconds (Continuous)
Atrial Tach	Atrial Tachycardia	160 BPM (Continuous)
Paroxysmal Atrial Tach	Paroxysmal Atrial Tachycardia	160 BPM for five seconds 80 BPM for ten seconds (Continuous)
Supravent Tach	Supraventricular Tachycardia	200 BPM (Continuous)
Sinus Arrhythmia	Sinus Arrhythmia	Normal beats at a fluctuating rate from 60 BPM to 100 BPM (Continuous)
Missed Beat – Auto	Missed Beat	NSR of 80 BPM with a missed beat (Missed Beat, 36 NSR) (Continuous)
Missed Beat – Man	Missed Beat	NSR of 80 BPM with a missed beat (One-Time Event)
Nodal Rhythm	Nodal Rhythm	60 BPM with very short P-R interval (Continuous)


VENTRICULAR		
Abbreviation	Arrhythmia	Description
Pair of PVCs – Auto	Pair of Premature Ventricular Contractions	NSR of 80 BPM with Periodic Group of 2 Type 1 PVCs (2 PVC Type 1, 36 NSR) (Continuous)
Pair of PVCs – Man	Pair of Premature Ventricular Contractions	NSR of 80 BPM with Periodic Group of 2 Type 1 PVCs (One-Time Event)
Run of 5 PVCs – Auto	Run of 5 Premature Ventricular Contractions	NSR of 80 BPM with periodic group of 5 Type 1 PVCs (5 PVC Type 1, 36 NSR) (Continuous)
Run of 5 PVCs – Man	Run of 5 Premature Ventricular Contractions	NSR of 80 BPM with periodic group of 5 Type 1 PVCs (One-Time event)
Run of 11 PVCs – Auto	Run of 11 Premature Ventricular Contractions	NSR of 80 BPM with periodic group of 11 Type 1 PVCs (11 PVC Type 1, 36 NSR) (Continuous)
Run of 11 PVCs – Man	Run of 11 Premature Ventricular Contractions	NSR of 80 BPM with periodic group of 11 Type 1 PVCs (One-Time event)
6 PVCs per Min	6 Premature Ventricular Contractions per minute	NSR of 80 BPM with 6 Type 1 PVCs per minute (Continuous)
12 PVCs per Min	12 Premature Ventricular Contractions per minute	NSR of 80 BPM with 12 Type 1 PVCs per minute (Continuous)
24 PVCs per Min	24 Premature Ventricular Contractions per minute	NSR of 80 BPM with 24 Type 1 PVCs per minute (Continuous)
Freq Multifocal PVCs	Frequent Multifocal Premature Ventricular Contractions	NSR of 80 BPM with every fourth beat being an alternating Type 1 and Type 2 PVC (Continuous)
Bigeminy	Bigeminal Rhythm	NSR of 80 BPM with every other beat a Type 1 PVC (Continuous)
Trigeminy	Trigeminal Rhythm	NSR of 80 BPM with every third beat a Type 1 PVC (Continuous)

Vent Tach	Ventricular Tachycardia	160 BPM, No P-wave, Beats similar to Type 1 PVC (Continuous)
Vent Fib – Coarse	Ventricular Fibrillation	Irregular waveform with no real P-wave or clear R-R interval and a high signal level (Continuous)
Vent Fib – Fine	Ventricular Fibrillation	Irregular waveform with no real P-wave or clear R-R interval and a low signal level (Continuous)
Asystole	Asystole	Flat line signal (Continuous)


CONDUCTION		
Abbreviation	Arrhythmia	Description
1st Deg Heart Block	First Degree Heart Block	80 BPM with a long P-R interval of 250 ms (Continuous)
2nd Deg Heart Block	Second Degree Heart Block	80 BPM with increasing P-R interval for four beats (160, 220, 400, 470 ms) followed by a P wave without a QRS (Continuous)
3rd Deg Heart Block	Third Degree Heart Block	80 BPM with P wave rate of 80 BPM and QRS rate of 30 BPM (Continuous)
Rt Bundle Branch Block	Right Bundle Branch Block	80 BPM with Normal P-wave and P-R interval but wider QRS complexes (Continuous)
Lf Bundle Branch Block	Left Bundle Branch Block	80 BPM with Normal P-wave and P-R interval but wider QRS complexes (Continuous)






ECG – PACEMAKER

The PS-2105 can send paced waveforms to ECG machines in 3, 5 or 12-lead configurations. It has independent outputs for each signal lead, referenced to the right leg.




There are 7 paced simulation signals available which model when a pacemaker accompanies the heartbeat. The PS-2105 is placed into PACEMAKER mode by pressing the  category key.

The display will resemble the following:



The pacemaker rhythms and signals can be selected by using   to highlight the parameter to change and using   to scroll to the desired option. Then  is used to accept the new setting.

Alternately, to see a submenu of all the options for a highlighted parameter, use .


Use   to scroll to the desired option. Then  is used to accept the new setting.

The following is a brief description of how the PS-2105 simulates the available Pacemaker Waveforms:


PACEMAKER		
Abbreviation	Waveform	Description
Atrial Pacer	Atrial Pacemaker Wave	80 BPM with Pacer Pulse at the start of each P wave
Asynchronous	Asynchronous Pacemaker Wave	75 BPM with Pacer Pulse at the start of each QRS wave and no P wave
Non-Capture	Ventricular Pacemaker Wave with Periodic Non-Response	75 BPM Ventricular Paced beats with every tenth beat not responding
Non-Function	Ventricular Pacemaker Wave with no Heart Response	75 BPM Ventricular Paced beats with no heart response
Demand – Occasional	Demand Pacemaker Wave with Occasional Sinus Beats	20 NSR beats followed by 20 Ventricular Paced beats
Demand – Frequent	Demand Pacemaker Wave with Frequent Sinus Beats	40 NSR beats followed by 40 Ventricular Paced beats
AV – Sequential	AV-Sequential Pacemaker Wave	75 BPM with Pacer Pulse at the start of both the P and QRS waves






ECG – PERFORMANCE





The PS-2105 can send performance waveforms to ECG machines in 3, 5 or 12-lead configurations. It has independent outputs for each signal lead, referenced to the right leg.

There are 15 Performance waves available for testing and verifying. The PS-2105 is placed into PERFORMANCE mode by pressing the  category key.

The display will resemble the following:



These waves and amplitudes can be selected by using   to highlight the parameter to change and using   to scroll to the desired option. Then  is used to accept the new setting.

Alternately, to see a submenu of all the options for a highlighted parameter, use . Use   to scroll to the desired option. Then  is used to accept the new setting.


Auto Wave

If the Performance parameter is set to AUTO, the unit will automatically sequence through all of the performance waves, starting with Square Wave .125 Hz, incrementing at a fixed interval. The interval may be set in the System Setup Menu under “Auto Step Time”.

A countdown timer is shown in the display:


Displays time (seconds) remaining before advancing to next waveform.

Performance	
Square Wave .125 Hz	
<hr/>	
Amplitude:	1.0 mV
<hr/>	
Auto Timer (28)	
<hr/>	


The  key can be used to exit the Auto Mode during the sequence.






TRAINING





The PS-2105 provides the unique feature of a Training Mode to aid the user in practicing the identification of arrhythmias. The unit will sequence through the arrhythmias, allowing the user to look at the output on their equipment, identify the arrhythmia and then verify their conclusion with the correct name shown on the display. The user can select either manual or timed sequencing, as well as whether the arrhythmias will display in order or randomly. Subsets of the Arrhythmias can be selected to allow for individualization.

The PS-2105 is placed into TRAINING mode by pressing the  category key.

The display will resemble the following:




The Timer, Randomizer and Arrhythmias can be selected by using   to highlight the parameter to change and using   to scroll to the desired option. Then  is used to accept the new setting.


Alternately, to see a submenu of all the options for a highlighted parameter, use . Use   to scroll to the desired option. Then  is used to accept the new setting.

When ready to begin the Training, use . The appropriate arrhythmia screen will be displayed with Training Mode indicated.

```
Premature Arrhythmias:
  Atrial PAC - Auto
ECG Artifact:      none
-----
Train:              All/Manual
(Enter→Next Quit→Pause)
```

ARTIFACT
None*
50 Hz
60 Hz
Muscle
Wander
(Baseline)
Resp
(Respiration)

If in the timed mode, the unit will switch to the next arrhythmia automatically at the set time. If in the manual mode, use  to go to the next arrhythmia when ready.

To pause or exit the Training Mode during a session, use . The following message box will be displayed:


```
Press ENTER to continue.
Press QUIT to exit.
```


SUBSET

The subset feature allows the user to select specific arrhythmias for a more controlled training. This feature is selected by setting the “Arrhythmias:” parameter to “Subset.”


The selection of the subset is done by marking those specific arrhythmias or groups of arrhythmias of interest. After a subset of arrhythmias has been selected, it will remain in memory. It may then be edited at any time prior to starting a training session.





The following procedure is used to modify the subset:

1. From the main screen, highlight Arrhythmias: and use  to open the choices screen.







2. Highlight Subset and use  to open the Arrhythmia Groups Submenu screen.




3. To select all the arrhythmias in a group, use   to scroll to the category and  to toggle the indicator to “ALL”. To select none of the arrhythmias in a group, use  to toggle the indicator to “NO”.

Arrhythmia Groups:	
NO	Premature
ALL	Supraventricular
ALL	Ventricular
ALL	Conduction
(ENTER → TOGGLE CHOICES → SUBSET)	


4. To select some of the arrhythmias in a group, use  to display the list of arrhythmias for a specific group. Then use   to scroll through the arrhythmias. Any arrhythmia marked with a check (✓) will be included in the subset. Use  to toggle the selection of an arrhythmia on and off.

Premature Arrhythmias:	
	Atrial PAC - Auto
	Nodal PNC - Auto
✓	PVC 1 - Auto
✓	PVC 1 Early - Auto
(Enter → Toggle Selection)	


NOTE: If the group is pre-selected with “ALL”, all of the arrhythmias will be checked, thus making it easy to deselect a few. If the group is pre-selected with “NO”, none of the arrhythmias will be checked, thus making it easy to select a few.

5. When completed selecting the desired arrhythmias from that group, use  to return to the Group Submenu. "SOME" will appear to indicate a partial selection of the arrhythmias in that group,

```
Arrhythmia Groups:
SOME Premature
NO  Supraventricular
ALL Ventricular
ALL Conduction
(ENTER → TOGGLE CHOICES → SUBSET)
```

Additional groups may be modified in the same manor. When done with all the groups, use  again to return to the Training Mode.






SETUP

The PS-2105 allows for setup of the System Parameters through the  category key.

The System Setup screen allows for the setting of the parameters controlling various function of the unit, as well as the viewing of Battery Life and Software information.

The display will resemble the following:

```
System Setup MORE↵
SpO2 Output      Disabled
Auto Off Timer (Min) 30
Backlight Time (Sec) 5
Backlight Intensity 100%
Battery Life      100%
```

These settings can be selected by using   to highlight the parameter to change and using   to scroll to the desired option. Then  is used to accept the new setting.

Use  or , , , , or  to exit from the setup screens.

The following is a brief description of the parameters and the available range of settings:

Parameter	Description	Range
SpO ₂ Output	Sets the ability to drive an external SpO ₂ module (MSP-2100)	Enabled/Disabled
Auto Off Timer	The elapsed time after which the unit will automatically power down. This timer is reset by each key depression. (Setting the value to 0 eliminates this feature.)	0-30 min
Backlight Timed	Off – Always off 1-30 sec – The elapsed time after which the backlight will automatically turn off. Manual – The backlight will be manually controlled by backlight key.	Off, 1-30 sec, Manual
Backlight Intensity	Sets the intensity of the backlight. (Note: Lower intensities extend battery life.)	0-100%
Battery Life	Displays current life of the batteries. At 5%, a warning screen will appear. At 0%, the unit will power down automatically.	5-100% (Read Only)
Contrast Adjust	Sets the contrast of the display screen.	0-20
Power up with	Selects the values that will be used when the unit is first turned on. It is also used to Set the Custom Defaults, if used. (See Power Up Settings).	Default/Last/Custom/ Set Custom Defaults
Auto Step Time	Sets the interval that is used with the Auto increment features in BPM and Performance.	1 to 60 sec
Software	Displays current software program.	(Read Only)

POWER UP SETTINGS

The PS-2105 allows the user to tailor the settings that the unit will have on Power Up. The “Power Up With” parameter in the System Setup Menu allows for the selection of either Default, Last or Custom selections.

Default

If this option is selected the following settings will be used every time the unit is turned on.

ECG – NSR: 80 BPM, 1.0 mV, Adult QRS, 0.0 mV ST Elevation, Artifact - None,

SpO₂ Output Disabled

ECG – Arrhythmia: 1.0 mV, Artifact - None, Premature - Atrial PAC - Auto

ECG – Performance: 2 Hz Square Wave, 1.0 mV

ECG – Pacemaker: Atrial, 5 mV Amplitude, 1.0 ms Width

SystemSetup:

Auto Timer Off	30 min
Backlight Time	5 sec
Backlight Intensity	100%
Contrast Adjust	10
Power Up With	Default
Auto Step Time	5 sec


Last

If this option is selected, the unit will remember the settings that were being used when it was turned off and bring them back when the power is turned on.

Custom

If this option is selected, the user may save a unique set of default parameters and the unit will recall them every time the power is turned on.

Set Current As Custom

To create the set of custom default parameters, this fourth choice is provided in this parameter. The user simply configures the unit to the desired default conditions, selects this option and presses . The current configuration is then saved as the Custom Power up values.

SpO₂ (Option)

The PS-2105 has the ability to drive an external SpO₂ module. This module (MSP-2100) accepts the FingerSim family of SpO₂ finger simulators (fingers are available with SpO₂ of 80, 90 and 97%). The output pulses the fingers at the NSR BPM rate (up to 180 BPM). The output is off in Arrhythmia and Performance Modes.

The module plugs directly into the AUX (7 pin mini din) connector and is powered from the PS-2105. The output is only functional when the unit is powered from the Battery Eliminator provided with the MSP-2100 Module, since the batteries do not have enough power to run this option.




The output is enabled and disabled in the System Setup screen.



OUTPUT CONNECTIONS

The following is a representation of the socket connector found on the unit. It is viewed as if looking at the socket in the unit, not the cable pins.

AUX CONNECTOR



MANUAL REVISIONS

<u>Revision #</u>	<u>Program #</u>	<u>Revisions Made</u>
Rev 01	DT7348	Origination
Rev 02	DT7348	Misc. Updates
Rev 03	DT7348	Overlay Additions
Rev 04	DT7348CC	Format Updated, Pictures Updated, Misc. Edits

LIMITED WARRANTY

WARRANTY: BC GROUP INTERNATIONAL, INC. WARRANTS ITS NEW PRODUCTS TO BE FREE FROM DEFECTS IN MATERIALS AND WORKMANSHIP UNDER THE SERVICE FOR WHICH THEY ARE INTENDED. THIS WARRANTY IS EFFECTIVE FOR TWELVE MONTHS FROM THE DATE OF SHIPMENT.

EXCLUSIONS: THIS WARRANTY IS **IN LIEU OF** ANY OTHER WARRANTY EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO ANY IMPLIED WARRANTY OF **MERCHANTABILITY** OR FITNESS FOR A PARTICULAR PURPOSE.

BC GROUP INTERNATIONAL, INC. IS NOT LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES.



NO PERSON OTHER THAN AN OFFICER IS AUTHORIZED TO GIVE ANY OTHER WARRANTY OR ASSUME ANY LIABILITY.

REMEDIES: THE PURCHASER'S SOLE AND EXCLUSIVE REMEDY SHALL BE: (1) THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS OR PRODUCTS, WITHOUT CHARGE. (2) AT THE OPTION OF **BC GROUP INTERNATIONAL, INC.**, THE REFUND OF THE PURCHASE PRICE.

SPECIFICATIONS

ECG SIMULATION			
RATE	NORMAL SINUS RHYTHM	30, 40, 45, 60, 80, 90, 100, 120, 140, 160, 180, 200, 220, 240, 260, 280, 300 BPM	
	PERFORMANCE WAVEFORMS	SINE	0.1, 0.5, 5, 10, 40, 50, 60, 100 Hz
		SQUARE	0.125, 2.0 Hz
		TRIANGLE	2.0, 2.5 Hz
	PULSE	30, 60, 120 BPM; 60 ms width	
ACCURACY	± 1%		
AMPLITUDE	0.5, 1.0, 1.5, 2.0 mV (Lead II)		
	ACCURACY	± 2% @ Lead II	
HIGH LEVEL	OUTPUT	200 times Amplitude	
	ACCURACY	± 5%	
QRS DURATION	ADULT	80 ms	
	PEDIATRIC	40 ms	
ST SEGMENT (ELEVATION)	± 0, 0.05, 0.1, 0.2, 0.3, 0.4, 0.5, 0.6, 0.7, 0.8 mV		
LEAD TO LEAD IMPEDANCE	1000 Ω		

PACEMAKER WAVEFORMS		
RATE	75 BPM	
	ACCURACY	± 1%
AMPLITUDE	1, 2, 3, 4, 5, 6, 7, 8, 9, 10 mV	
	ACCURACY	± 10%
WIDTH	0.1, 0.5, 1.0, 1.5, 2.0 ms	
	ACCURACY	± 5%

PHYSICAL & ENVIRONMENTAL		
DISPLAY	LCD Graphical 128 X 64 Pixels, White LED Backlight	
CONSTRUCTION	ENCLOSURE	ABS Plastic
	FACE PLATE	Lexan, Back printed
SIZE	8.80 x 6.04 x 1.72 Inches (223.5 x 153.4 x 43.7 mm)	
WEIGHT	< 2 Lbs (0.91 kg)	
OPERATING RANGE	15 to 40 °C (59 to 104 °F)	
STORAGE RANGE	-20 to 65 °C (-4 to 149 °F)	
ELECTRICAL		
BATTERY	9V Alkaline Battery (2 Required) (ANSI/NEDA 1604A or equivalent)	
BATTERY ELIMINATOR	WITHOUT MSP-2100	9 VDC, 200 mA  BC20-21100 (USA Version) BC20-21101 (Euro Version)
	WITH MSP-2100	10 VDC, 500 mA  BC20-21103 (USA Version) BC20-21101 (Euro Version)

NOTES



**BC GROUP INTERNATIONAL, INC.
3081 ELM POINT INDUSTRIAL DRIVE
ST. CHARLES, MO 63301
USA**

1-800-242-8428

1-314-638-3800

www.bcgrouptl.com
sales@bcgrouptl.com

**PS-2105 Patient Simulator User Manual
09/12 – Rev 04**

**Copyright © 2012
Made in the USA**